



Venturi valve Intersurgical EcoLite™ mask kits



Oxygen and Aerosol Therapy ▪ Fixed Oxygen Concentration (High Flow)



Quality, innovation and choice

Interact with us



www.intersurgicalinc.com



Intersurgical EcoLite™ venturi valves and mask kits

Introducing a comprehensive range of venturi valves available individually or in kits. The venturi valves are designed to meet the requirements of BS EN 13544-3:2001+A1:2009 and are color coded for easy identification.

Optimum flow rates

The flow printed on the valves ensures that in most cases, the patient receives a total flow of approximately 40 L/min. This is designed to meet the patient's 'normal' peak flow requirements.

In situations where a patient's peak flow requirements may exceed 40 L/min, the oxygen flow can be increased leading to a higher total flow (see table) with no significant change in delivered oxygen concentration.

Flow rates												
Oxygen flow meter — Liters per minute (L/min)												
Oxygen Concentration	2	3	4	5	6	7	8	9	10	11	12	15
	24%	52	78	104								
	28%			45	56	67	79	90				
	31%					47	55	63	71	79	87	
	35%						45	51	56	62	68	84
	40%								42	46	50	62
60%												30
Total flow rate to patient — Liters per minute (L/min)												

NOTE: Figures in **bold** indicate total flow rate to patient when using recommended oxygen flow.

References:

- 1. Clinical Practice in Respiratory Care, James B. Fink & Gerald E. Hunt, 1999
- 2. British Thoracic Society Guideline for oxygen use in adults in healthcare and emergency settings, O'Driscoll BR, et al. Thorax 2017
- 3. IQR 226993
- 4. IQR 319281



Welcome to the Comfort Zone

The unique design and materials used to manufacture the Intersurgical EcoLite™ mask provides improved patient comfort and a reduced environmental impact.

Comfortable for the Patient

Improved patient comfort has been key to the development of the Intersurgical EcoLite mask. The latest manufacturing technology has enabled us to combine two non-PVC materials in the same product, the polypropylene material forming the body of the Intersurgical EcoLite mask is clear, lightweight and rigid enough to maintain the masks shape, while the second softer TPE material is utilized in the seal, which is in contact with the patient's face.

Comfortable for the Environment

The Intersurgical EcoLite mask is an important part of our Eco range, designed as part of our ongoing focus on sustainable development. Various groups have questioned the use of PVC in medical products and its impact on the environment. The utilization of alternative materials has resulted in the elimination of PVC from the Intersurgical EcoLite mask body and seal, reducing the environmental impact.

Intersurgical EcoLite

[Find out more](#)

Eco range

[Find out more](#)



Incurved nose seal

Conforms to different nose shapes designed to reduce the risk of oxygen entering patient's eyes

No metal nose clip

MRI compatible

Soft face seals

Increased patient comfort


Elastic can be positioned under or over the ears


Below ear position eliminates trauma to top of ears



A choice of "under chin" or "on chin" positions

Provides a better fit on a wider range of patient face shapes





Code	Description	Color	Tube length	Box Qty.
1060085	Intersurgical EcoLite™, adult, oxygen mask with 60% venturi valve, green			40
1060081	Intersurgical EcoLite, adult, oxygen mask with 60% venturi valve, green and tube		1.8 m	40
1040085	Intersurgical EcoLite, adult, oxygen mask with 40% venturi valve, red			40
1040081	Intersurgical EcoLite, adult, oxygen mask with 40% venturi valve, red and tube		1.8 m	40
1035085	Intersurgical EcoLite, adult, oxygen mask with 35% venturi valve, yellow			40
1035081	Intersurgical EcoLite, adult, oxygen mask with 35% venturi valve, yellow and tube		1.8 m	40
1031085	Intersurgical EcoLite, adult, oxygen mask with 31% venturi valve, orange			40
1031081	Intersurgical EcoLite, adult, oxygen mask with 31% venturi valve, orange and tube		1.8 m	40
1028085	Intersurgical EcoLite, adult, oxygen mask with 28% venturi valve, white			40
1028081	Intersurgical EcoLite, adult, oxygen mask with 28% venturi valve, white and tube		1.8 m	40
1024085	Intersurgical EcoLite, adult, oxygen mask with 24% venturi valve, blue			40
1024081	Intersurgical EcoLite, adult, oxygen mask with 24% venturi valve, blue and tube		1.8 m	40
1107085	Intersurgical EcoLite, adult, venturi mask kit		1.8 m	25
1107012	Intersurgical EcoLite, pediatric, venturi mask kit		1.8 m	30

Code	Description	Color	Box Qty.	Code	Description	Color	Box Qty.
1400060	Venturi valve 60% oxygen, green		50	1400031	Venturi valve 31% oxygen, orange		50
1400050	Venturi valve 50% oxygen, pink		50	1400028	Venturi valve 28% oxygen, white		50
1400040	Venturi valve 40% oxygen, red		50	1400024	Venturi valve 24% oxygen, blue		50
1400035	Venturi valve 35% oxygen, yellow		50	1040013	Oxygen recovery kit, 40% oxygen		40

Make an inquiry

Lower environmental impact product

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 298.2 / IS10.22 US • Issue 16 12.24



Intersurgical Inc, 2055 Executive Dr. Indianapolis, IN 46241, USA

T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • Ireland • Italy • Japan
Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA