



Universal Connector Oxygen Tubing

Intersurgical has combined the Universal Connector with our Intersurgical EcoLite™ Oxygen and Aerosol Therapy product portfolio



Oxygen and Aerosol Therapy ▪ Oxygen Therapy



Universal Connector oxygen tubing

Intersurgical has combined the Universal Connector oxygen tubing with our Intersurgical EcoLite Oxygen and Aerosol Therapy product portfolio.

This combination of the two product lines can improve efficiency for the clinician and provide added safety for the patient. The Universal Connector allows the Intersurgical EcoLite Product Range to be connected to the flow meter without the use of a traditional "Christmas Tree" adapter. This innovative combination helps eliminate the need to search for or replace oxygen flow meter "Christmas Tree" adapters.

Patient Safety

Helping to eliminate the safety concern of not being able to connect to flow meters because of missing adapters.



Comfortable for the environment

Various groups have questioned the use of PVC in medical products and its impact on the environment. Intersurgical's aim is to reduce the environmental impact of its products and processes. The utilization of new materials and the latest manufacturing technology has resulted in the elimination of PVC from the mask shell resulting in a reduced environmental impact.

Incurved nose seal

Conforms to different nose shapes designed to prevent oxygen entering patient's eyes

No metal nose clip

MRI compatible

Exhalation ports

Designed to minimize the amount of rebreathed carbon dioxide

A choice of "under chin" or "on chin" positions

Provides a better fit on a wider range of patient face shapes

Crush resistant tubing


Allows oxygen to continue to flow

Provides a direct connection

Eliminates the need for an additional adapter or "Christmas tree"



us.intersurgical.com/info/intersurgicalecolite



Lighter

50% lighter than a PVC mask, more comfortable for the patient

Greener


The elimination of PVC from the masks has dramatically lowered their environmental impact

Soft face seals

Increased patient comfort

Elastic can be positioned under or over the ears

Below ear position eliminates trauma to top of ears



Code	Description	Tube length	Box Qty.
11366100	Intersurgical EcoLite, adult, medium concentration oxygen mask with Universal Connector tube	1.8m	30
11376100	Silente™ Eco, adult, adjustable venturi mask kit with Universal Connector tube	1.8m	50
11386100	Silente, Intersurgical EcoLite, pediatric adjustable venturi mask kit with Universal Connector tube	1.8m	30
11426100	Sentri™ Intersurgical EcoLite, adult, mask with CO ₂ monitoring line with Universal Connector tube	1.8m	30
11816100	Intersurgical EcoLite, adult, high concentration oxygen mask with Universal Connector tube	1.8m	20
11906100	Intersurgical EcoLite, pediatric, medium concentration oxygen mask with Universal Connector tube	1.8m	40

Lower environmental impact product

IS10.16 US • Issue 8 04.18