



Economy anesthetic face masks



Anesthesia ▪ Anesthetic Face Masks



Economy anesthetic face masks

A range of four color-coded economy PVC anesthetic face masks with a soft cushion and clear housing. Four sizes from small adult to extra large adult.



Clear flexible housing and cushion

For good patient visibility

Soft cushion

For patient comfort and air tight fit







Color-coded hook ring

For ease of identification, and integration with other Intersurgical airway products

Four sizes available

From small adult to extra large adult

Code	Description	Color	Box Qty.
1517000	Economy, anesthetic face mask, size 6, extra large adult, 22ID		25
1516000	Economy, anesthetic face mask, size 5, large adult, 22ID		30
1515000	Economy, anesthetic face mask, size 4, medium adult, 22ID		40
1514000	Economy, anesthetic face mask, size 3, small adult, 22ID		50

Make an inquiry