



i-flo™ high flow nasal cannula



Critical Care ▪ i-flo high flow nasal cannula



i-flo™ high flow nasal cannula


The i-flo is Intersurgical's single use patient interface designed for the delivery of High Flow Oxygen Therapy (HFOT). The i-flo is available in three sizes and is for use with adult patients in a hospital environment.

HFOT has become a key mode of respiratory support in patients with hypoxemic respiratory failure and is considered to offer a number of benefits versus more conventional oxygen therapies. These benefits have been shown to include:

- Improved oxygenation
- Reduced work of breathing
- Improved secretion clearance
- Improved patient tolerance
- Avoidance of intubation^[1]

The i-flo complements Intersurgical's extensive range of Critical Care breathing systems and interfaces and the unobtrusive design has been developed with a focus on performance and patient comfort.


The i-flo includes a number of unique features and benefits to help ensure the optimal delivery of high flow oxygen therapy.



[1] Masaji Nishimura. High-Flow Nasal Cannula Oxygen Therapy in Adults: Physiological Benefits, Indication, Clinical Benefits, and Adverse Effects. Resp.Care 2016 Apr;61(4):529-41



Adaptive and secure



The i-flo's unique, 180° swivel elbow enables easy and convenient positioning of bedside equipment.



The i-flo clip and lanyard tube management system securely supports the breathing system tubing and allows the patient to move freely without concerns of cannula dislocation or discomfort.

Comfortable and supportive



An anatomical face-mount design directs the head strap above the patient's ears, helping to minimize rubbing and discomfort.

The flexible, semi-rigid face-mount also ensures that strap pressure is distributed evenly over the patient's cheeks. This increases cannula stability and support and helps to prevent pressure hot-spots from developing.





i-flo™ high flow nasal cannula

Code	Description	Box Qty.
7300000	i-flo high flow nasal cannula, small adult	7
7301000	i-flo high flow nasal cannula, medium adult	7
7302000	i-flo high flow nasal cannula, large adult	7


Make an inquiry

Watch the video



Single limb breathing systems

Code	Description	Box Qty.
2026035	22 mm Flextube single limb heated wire with limb, ≥ 1.6 m length	20
2026036	22 mm Flextube single limb heated system with auto-fill chamber and limb, ≥ 1.6 m length	20



Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 251.2 / IS1.53 US ▪ Issue 4 01.25



Intersurgical Inc, 2055 Executive Drive, Indianapolis, IN 46241, USA
T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK ▪ Australia ▪ Benelux ▪ Brazil ▪ Canada ▪ China ▪ Colombia ▪ Czech Republic ▪ Denmark ▪ France ▪ Germany ▪ Ireland ▪ Italy ▪ Japan
Lithuania ▪ Malaysia ▪ Mexico ▪ Norway ▪ Philippines ▪ Portugal ▪ Russia ▪ South Africa ▪ Spain ▪ Sweden ▪ Taiwan ▪ Turkey ▪ USA