



Right-angled water traps

Self-sealing water traps



Critical Care • Accessories



Right-angled water traps

Self-sealing water traps

Introducing Intersurgical's new right-angled water traps for the collection of condensation.

Complementing our existing range of products, the water trap is now available in two sizes. Medium volume 70ml and high volume 125ml are both ideal for use in applications where a 90° attachment may be required.

The self-sealing facility allows the water trap to be emptied without interrupting flow.

Both water traps feature a 22ID machine-side connection to a 22OD breathing system-side connection, and can be attached directly to a ventilator or to an anesthetic machine.



1916000

1915000

Ordering information

| Code | Description | Box Qty. |
|---------|--|----------|
| 1916000 | Water trap, right angled, medium 70ml | 55 |
| 1915000 | Water trap, right angled, large 125ml | 40 |