



Right-angled water traps

Self-sealing water traps



Critical Care ▪ Accessories



Quality, innovation and choice

Interact with us



www.intersurgicalinc.com



Right-angled water traps

Self-sealing water traps

Intersurgical's right-angled water traps for the collection of condensation compliments our existing range of products, the angled water trap is available in two sizes.

Medium volume and large volume, both are ideal for use in applications where a 90° attachments may be required.

The self-sealing facility allow the water traps to be emptied without interrupting gas flow.

Both water traps feature a 22ID machine side connection to a 22OD breathing circuit side connection, and can be attached directly to a ventilator or to an anesthetic machine.



Self-sealing facility

Allow the water traps to be emptied without interrupting flow



Code	Description	Box Qty.
1916000	Water trap, self-sealing, right-angled, medium volume	55
1915000	Water trap, self-sealing, right-angled, large volume	40

Make an enquiry

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 131.2 / IS1.33 US • Issue 6 10.24



Intersurgical Inc, 2055 Executive Dr. Indianapolis, IN 46241, USA

T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • Ireland • Italy • Japan
Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA