

AccuPAP®

Positive Airway Pressure Therapy System







Oxygen and Aerosol Therapy ▪ Pulmodyne®

INTERSURGICAL®
COMPLETE RESPIRATORY SYSTEMS

Quality, innovation and trust


Interact with us

www.intersurgicalinc.com



AccuPAP®

The AccuPAP system is designed to deliver positive airway pressure indicated for the prevention or reversal of atelectasis in adult patients. The AccuPAP includes adjustable PEEP, this allows the healthcare professional to deliver pressure levels according to each patient's requirements. The PEEP can be adjusted without requiring a change to be made to the input gas flow.



Code	Description	Box Qty.
313-6000	AccuPAP with mouthpiece, standard tubing	10
313-6001	AccuPAP with mouthpiece and manometer, standard tubing	10
313-6200	Apep PEP therapy system	10
313-6002	AccuPAP with mouthpiece, universal tubing	10
313-6003	AccuPAP with mouthpiece and manometer, universal tubing	10

Make an inquiry Watch the video

AccuPAP and DPM are manufactured by Pulmodyne Inc, and distributed by Intersurgical Ltd.

U.S. Lit No. 379.1 / IS10.40 US ▪ Issue 7.1 05.26

INTERSURGICAL® COMPLETE RESPIRATORY SYSTEMS

Intersurgical Inc, 2055 Executive Drive, Indianapolis, IN 46241, USA
T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com